JenBiOxzyme®

X S - E X P L O R E R - T E S T K I T

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Puffer: Puffer 1: Mangan-Peroxidase, Versatile Peroxidase

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Activity: 12.6 U/mg

26.3 U/mg protein

Unit definition: One unit will oxidize 1.0 μmole of Mn(II) per minute at pH 4.5 at 25°C.

Assay method: The reaction velocity of the enzyme-catalyzed reaction is determined by the increase in absorbance at 270 nm.

Form: Lyophilisate, containing 477 μg protein/mg (BCA-Assay, Pierce)

Appearance: Brownish powder

Impurities: Laccase: 0.16 U/mg

Lignin peroxidase: 0.09 U/mg

Stability: Stable up to 12 months, when stored at -18°C

pH-Optimum: Mn(II) oxidation (50 mM malonate buffer): pH 4.5

Notes: This hemoprotein is involved in the oxidative degradation of lignin in white-rot basidiomycetes.

Product No.: EN- 202

Lignin-Peroxidase from *Phanerochaete chrysosporium*

Diarylpropane:oxygen, hydrogen-peroxide oxidoreductase (EC 1.11.1.14 , CAS 93792-13-3)

Lot No.: P-16-JB03

Source: *Phanerochaete chrysosporium*

Reaction: $3,4\text{-dimethoxybenzyl alcohol} + \text{H}_2\text{O}_2 \longrightarrow 3,4\text{-dimethoxybenzaldehyde} + 2 \text{H}_2\text{O}$

Activity: 12.6 U/mg

26.3 U/mg protein

Unit definition: One unit will oxidize 1.0 μmole of 3,4-dimethoxybenzyl alcohol (veratryl alcohol) per minute at pH 3.3 at 25°C

Assay method: The reaction velocity of the enzyme-catalyzed reaction is determined by the increase in absorbance at 310 nm.

Form: Freeze-dried powder containing 298 μg protein/mg (BCA-Assay, Pierce)

Appearance: Brownish powder

Impurities: Laccase: 0.16 U/mg

Manganese Peroxidase: 0.09 U/mg

Stability: Stable up to 12 months, when stored at -18°C

pH-Optimum: 3,4-dimethoxybenzyl alcohol oxidation (100 mM sodium tartrate buffer): pH 3.3

Notes: A hemoprotein, involved in the oxidative breakdown of lignin by white-rot basidiomycete fungi. The enzyme from *P. chrysosporium* oxidizes typical peroxidase dye substrates at the heme iron, a reaction involving the formation of compound II (FeIV=O); it also oxidizes 3,4-dimethoxybenzyl alcohol (veratryl alcohol) to the radical cation. The bound veratryl alcohol radical is proposed to bring about the oxidative cleavage of C-C and ether (C-O-C) bonds in lignin model compounds of the diarylpropane and arylpropane-aryl ether type.

Product No.: EN- 203

Versatile Peroxidase from *Bjerkandera adusta*

Mn(II): hydrogen-peroxide oxidoreductase (EC 1.11.1.16, CAS 42613-30-9, 114995-15-2)

Lot No.: U-09-JB53

Source: *Bjerkandera adusta*

Reaction: $2 \text{ Mn(II)} + 2 \text{ H}^+ + \text{ H}_2\text{O}_2 \longrightarrow 2 \text{ Mn(III)} + 2 \text{ H}_2\text{O}$
 $3,4\text{-dimethoxybenzyl alcohol} + \text{ H}_2\text{O}_2 \longrightarrow 3,4\text{-dimethoxybenzyl aldehyde} + 2 \text{ H}_2\text{O}$

Activity: 3.6 U/mg Lyo. Mn(II)-oxidizing activity; 15.6 U/mg protein
 0.2 U/mg Lyo. 3,4-dimethoxybenzyl alcohol-oxidizing activity; 0.8 U/mg protein

Unit definition: The amount of one unit will oxidize 1.0 μmole of Mn(II) per minute at pH 4.5 at 25°C, or 0.05 μM of 3,4-dimethoxybenzyl alcohol per minute at pH 3.0 (25°C).

Assay method: The reaction velocity of the enzyme-catalyzed reaction is determined by the increase in absorbance at 270 nm resulting from the formation of Mn(III)-malonate, or in absorbance at 310 nm resulting from the formation of 3,4-dimethoxybenzyl aldehyde.

Form: Freeze-dried powder, containing 232 μg protein/mg (BCA-Assay, Pierce)

Appearance: Brownish powder

Impurities: Lignin peroxidase: 0.19 U/mg

Stability: Stable up to 12 months, when stored at -18°C

pH-Optimum: Mn(II) oxidation (50 mM malonate buffer): pH 4.5; 3,4-dimethoxybenzyl alcohol (veratryl alcohol) (100 mM tartrate buffer): pH 3.0

Notes: Versatile peroxidase (syn. manganese-lignin peroxidase) is a new ligninolytic enzyme, combining catalytic properties of manganese peroxidase (oxidation of Mn(II)), lignin peroxidase (Mn-independent oxidation of non-phenolic aromatic compounds) and plant peroxidase (oxidation of hydroquinones and substituted phenols).

Product No.: EN- 204

Laccase from *Cerrena unicolor*

Benzenediol:oxygen oxidoreductase (EC 1.10.3.2, CAS 80498-15-3)

Lot No.: C-16-JB04

Source: *Cerrena unicolor*

Reaction: $4 \text{ benzenediol} + \text{ O}_2 \longrightarrow 4 \text{ benzosemiquinone} + 2 \text{ H}_2\text{O}$

Activity: 123 U/mg
 123 U/mg protein

Unit definition: One unit will oxidize 1.0 μmole of ABTS per minute at pH 4.5 and 25°C.

Assay method: The velocity of the enzyme-catalyzed reaction is determined by the increase in absorbance at 420 nm.

Form: Freeze-dried powder, containing 1234 μg protein/mg (BCA-Assay, Pierce)

Appearance: Brownish powder

Impurities: Mn-Peroxidase: 123 U/mg

Stability: Stable up to 12 months, when stored at -18°C

pH-Optimum: ABTS-oxidation (50 mM citrate/phosphate buffer): pH 4.5

Notes: A group of multi-copper proteins of low specificity acting on both o- and p-quinols, and often acting also on aminophenols and phenylenediamine. The semiquinone may react further either enzymically or non-enzymically.

Literatur

Alberts JF, Gelderblom WCA, Botha A, van Zyl WH (2009): Degradation of aflatoxin B1 by fungal laccase enzymes. *International Journal of Food Microbiology* 135: 47-52.

Arora DS, Sharma RK (2010): Ligninolytic Fungal Laccases and Their Biotechnological Applications. *Applied biochemistry and biotechnology* 160 (6): 1760-1788

Conesa A, Punt PJ, van den Hondel C (2002): Fungal peroxidases: molecular aspects and applications. *Journal of biotechnology* 93 (2): 143-158.

Franssen MCR, Steunenbergh P, Scott EL, Zuilhof H, Sanders JPM (2013): Immobilised enzymes in biorenewables production. *Chemical Society Reviews* 42 (15): 6491-6533.

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Gassara F, Brar SK, Verma M, Tyagi RD (2013): Bisphenol A degradation in water by ligninolytic enzymes. *Chemosphere* 92 (10): 1356-1360.

Fernando T, Bumpus JA, Austi SD (1990): Biodegradation of TNT (2,4,6-Trinitrotoluene) by *Phanerochaete chrysosporium*. *Applied and Environmental Microbiology* 1666-1671.

Hatakka (2004) Biodegradation of Lignin. (in Steibüchel [Ed.] *Biopolymers*. Wiley-VHC).

Hadibarata T, The ZC, Rubiyatno, Zubir MMFA, Khudhair AB, Yusoff ARM, Salim MR, Hidayat T (2013): Identification of naphthalene metabolism by white rot fungus *Pleurotus eryngii*. *Bioprocess Biosyst. Eng.* 36: 1455-1461.

Hadibarata T, Kristani RA (2014): Potential of a white-rot fungus *Pleurotus eryngii* F032 for degradation and transformation of fluorine. *Fungal Biology* 118: 222-227.

Hofrichter M, Wesenberg D, Rogalski J (2001): Fungal laccase: properties and activity on lignin. *J. Basic Microbiol.* 41 (3-4): 185-227.

Hofrichter M (2002): Review: lignin conversion by manganese peroxidase (MnP). *Enzyme and microbial technology* 30 (4): 454-466.

Hofrichter M, Ullrich R (2010a): *New Trends in Fungal Biooxidation*. In *Industrial Applications*, Hofrichter M, Ed. Springer Berlin Heidelberg: Vol. 10, pp 425-449.

Hofrichter M, Ullrich R, Pecyna MJ, Liers C, Lundell T (2010b): New and classic families of secreted fungal heme peroxidases. *Applied microbiology and biotechnology* 87 (3): 871-97.

Knezevic A, Milovanovic I, Stajic M, Loncar N, Brceski I, Vukojevic J, Cilerdzic J (2013): Lignin degradation by selected fungal species. *Bioresource Technology* 138: 117-123.

Knop D, Ben-Ari J, Salame TM, Levinson D, Yarden O, Hadar Y (2014): Mn²⁺-deficiency reveals a key role for the *Pleurotus ostreatus* versatile peroxidase (VP4) in oxidation of aromatic compounds. *Appl Microbiol Biotechnol* 98: 6795-6804.

Koschorreck K, Richter SM, Swierczek A, Beifuss U, Schmid RD, Urlacher VB (2008): Comparative characterization of four laccases from *Trametes versicolor* concerning phenolic C-C coupling and oxidation of PAHs. *Archives of Biochemistry and Biophysics* 474: 213-219.

Kullmann SW, Matsumura F (1996): Metabolic Pathways Utilized by *Phanerochaete chrysosporium* for Degradation of the Cyclodiene Pesticide Endosulfan. *Biology* 593-600.

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Lloret L, Eibes G, Feijoo G, Moreira MT, Lema JM (2012): Degradation of estrogens by laccase from *Myceliophthora thermophila* in fed-batch and enzymatic membrane reactors. *Journal of hazardous materials* 213-214: 175-183.

Lundell TK, Makela MR, Hilden K (2010): Lignin-modifying enzymes in filamentous basidiomycetes - ecological, functional and phylogenetic review. *Journal of basic microbiology* 50 (1): 5-20.

Martinez AT, Speranza M, Ruiz-Duenas FJ, Ferreira P, Camarero S, Guillen F, Martinez MJ, Gutierrez A, del Rio JC (2005): Biodegradation of lignocellulosics: microbial chemical, and enzymatic aspects of the fungal attack of lignin. *International Microbiology* 8 (3): 195-204.

Mayer AM, Staples RC (2002): Laccase: new functions for an old enzyme. *Phytochemistry* 60 (6): 551-565.

Minussi RC, de Moraes SG, Pastore GM, Duran N (2001): Biodecolorization screening of synthetic dyes by four white-rot fungi in a solid medium: possible role of siderophores. *Letters in Applied Microbiology* 33: 21-25.

Nakayama T, Amachi T (1999): Fungal peroxidase: its structure, function, and application. *Journal of Molecular Catalysis B-Enzymatic* 6 (3): 185-198.

Novotný C, Erbanová P, Sasek V, Kubátová A, Cajthaml T, Lang E, Krahl J, Zadrazil F. (1999): Extracellular oxidative enzyme production and PAH removal in soil by exploratory mycelium of white rot fungi. *Biodegradation* 10: 159-168.

Okeke C, Paterson A, Smith JE, Watson-Craik IA (1997): Comparative biotransformation of pentachlorophenol in soils by solid substrate cultures of *Lentinula edodes*. *Appl. Microbiol. Biotechnol.* 48: 563-569.

Pointing B (2001): Feasibility of bioremediation by white-rot fungi. *Appl. Microbiol. Biotechnol.* 57: 20-33.

Ruiz-Duen FJ, Morales M, Garcia E, Miki Y, Martinez MJ, Martinez AT (2009): Substrate oxidation sites in versatile peroxidase and other basidiomycete peroxidases. *Journal of Experimental Botany* 60: 441-452.

Salvachua D, Prieto D, Mattinen ML, Tamminen T, Liitia T, Lille M, Willfor S, Martinez AT, Martinez MJ, Faulds CB (2013): Versatile peroxidase as a valuable tool for generating new biomolecules by homogeneous and heterogeneous cross-linking. *Enzyme and Microbial Technology* 52: 303-311.

Schmid RD, Urlacher V (2007): *Modern Biooxidation*. Wiley.

da Silva Coelho-Moreira J, Bracht A, da Silva de Souza AC, Oliveira RF, de Sá-Nakanishi AB, Marques de Souza CG, Peralta RM (2013): Degradation of Diuron by *Phanerochaete chrysosporium*: Role of Ligninolytic Enzymes and Cytochrome P450. *BioMed Research International* ID 251354. 9 p.

Torres-Duarte C, Roman R, Tinoco R, Vazquez-Duhalt R (2009): Halogenated pesticide transformation by a laccase–mediator system. *Chemosphere* 77: 687-692.

Tsujiyama S, Muraoka T, Takada N (2013): Biodegradation of 2,4-dichlorophenol by shiitake mushroom (*Lentinula edodes*). *Biotechnol. Lett.* 35: 1079-1083.

Vasquez C, Anderson D, Oyarzun M, Carvajal A, Palma C (2014): Method for the stabilization and immobilization of enzymatic extracts and its application to the decolorization of textile dyes. *Biotechnol. Lett.* 36 (10): 1999-2010.

Wen X, Jia Y, Li J (2010): Enzymatic degradation of tetracycline and oxytetracycline by crude manganese peroxidase prepared from *Phanerochaete chrysosporium*. *Journal of Hazardous Materials* 177: 924-928.

Wen X, Jia Y, Li j (2009): Degradation of tetracycline and oxytetracycline by crude lignin peroxidase prepared from *Phanerochaete chrysosporium* – A white rot fungus. *Chemosphere* 75: 1003-1007.

Wong DWS (2009): Structure and Action Mechanism of Ligninolytic Enzymes. *Appl. Biochem. Biotechnol.* 157: 174-209.

Yehia RS (2014): Aflatoxin detoxification by manganese peroxidase purified from *Pleurotus ostreatus*. *Brazilian Journal of Microbiology* 45 (1): 127-133.

Zacchi L, Cox JR, Cheke RA, van der Walt E, Harvey PJ (1999): Biodegradation of Fenthion by *Phanerochaete Chrysosporium*. Workshop on Research Priorities for Migrant Pests of Agriculture in Southern Africa, Plant Protection Research Institute, Pretoria, South Africa, 24–26 March 1999.

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Hydrolasen (u.a. Esterase, Glycosidase, Amylase)*

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Single Fungal Enzymes

The JenaBios GmbH is able to produce in a position together with its partners a large number of enzymes from the classes of oxidoreductases and hydrolases on wild strain fermentation (see below). This extracellular multifunctional biocatalysts are able to carry out diverse biochemical changes in complex structures such as lignin or xenobiotics (e.g., Antibiotics, pesticides or drugs). They catalyze the oxidation of non-phenolic lignin and phenolic compounds. These fungal enzymes needs for the biotransformation reactions some co-substrates, who are economically priced and very simple to use (e.g. hydroperoxides such H₂O₂) or manganese).

Other key points of our product portfolio consisting the analysis, characterization of enzyme activities on reference substrates, as well as scale-up processes of new fermentation processes (up to 30-L scale).

We are produce our fungal enzymes from different organisms as freeze dried powder. The standard packages are 200 and 1,000 units (bulk amounts on request).

Available Enzymes	Cat.No.
Manganese-Peroxidase (<i>Nematoloma frowardii</i>)	EN-201
Lignin-Peroxidase (<i>Phanerochaete chrysosporium</i>)	EN-202
Versatile Peroxidase (<i>Bjerkandera adusta</i>)	EN-203
Laccase (<i>Cerrena unicolor</i> , <i>Trametes versicolor</i> *)	EN-204

Further Enzymes

Peroxygenases (*Agrocybe aegerita*, *Marasmius rotula*)*

Alcohol reductases*

Hydrolases (e.g.. Esterase, Glycosidase, Amylase)*

*) will be produced on demand

Assays

Enzymatic Assay of Manganese- and Versatile Peroxidase

Principle: Manganese peroxidase catalyzes the oxidation of Mn(II) to Mn(III) by H₂O₂. The highly reactive Mn(III) is stabilized via chelation in the presence of dicarboxylic acid [1].

Versatile peroxidase catalyzes the oxidation of Mn(II) to Mn(III) by H₂O₂. The highly reactive Mn(III) is stabilized via chelation in the presence of dicarboxylic acid [1].

Reaction: $2 \text{Mn(II)} + 2 \text{H}^+ + \text{H}_2\text{O}_2 \longrightarrow 2 \text{Mn(III)} + 2 \text{H}_2\text{O}$

Method: Continuous spectrophotometric determination of the increase in absorbance at 270 nm resulting from the formation of Mn(III)-malonate ($\epsilon_{270\text{nm}} = 11,590 \text{ M}^{-1} \text{ cm}^{-1}$) [2].

Conditions: Temperature: 25°C
pH: 4.5
Wavelength: 270 nm
Light path: 1 cm

Unit definition: The amount of enzyme that oxidizes 1.0 µmole of Mn(II) per minute at pH 4.5 and 25°C.

Reagents: Sodium-malonate buffer 50 mM, pH 4.5
H₂O₂ 20 mM
MnCl₂ 20 mM
Enzyme solution ca. 0.5 U/ml

All reagents should be prepared from deionized, distilled water. Dissolve the enzyme at a concentration of 0.3 mg/ml in distilled water. Dilute further in distilled water if necessary.

Procedure: *Pipette into a suitable quartz cuvette:*

buffer	690 µl
MnCl ₂ -solution	100 µl
Enzyme solution	200 µl

Incubate at 25°C and check the temperature, place the reaction mixture into the photometer and monitor the absorbance until constant, start the enzymatic reaction by adding

H ₂ O ₂	10 µl
-------------------------------	-------

Immediately mix and record the increase in absorbance at 270 nm for ca. 1 minute. Calculate the $\Delta A/\text{min}$ from the linear range of the curve (initial range: standardized as 10 sec after start).


Calculation:

$$\begin{aligned} \text{Volume activity} &= \frac{\Delta A/\text{min} \cdot 1000 \cdot \text{volume}_{\text{reaction mixture}} (\mu\text{l})}{11,590 \cdot \text{volume}_{\text{enzyme solution}} (\mu\text{l})} \\ &= [\text{U/ml enzyme solution}] \end{aligned}$$

References: [1] Hofrichter M *et al.* (1998). *Appl. Environ. Microbiol.* 64: 399-404.
[2] Wariishi H *et al.* (1992). *J. Biol. Chem.* 267: 23688-23695.

Enzymatic Assay of Laccase

Principle: Fungal laccase oxidizes 2,2'-azino-bis(3-ethylbenzthiazoline-6-sulfonate) (ABTS) to the stable green-colored cation radical (ABTS^{•+}).

Method: Continuous spectrophotometric determination of the increase in absorbance at 420 nm [1;2].

Conditions: Temperature: 25°C
 pH: 4.5
 Wavelength: 420 nm
 Light path: 1 cm

Unit definition: The amount of enzyme that oxidizes 1.0 μmole of ABTS per minute at pH 4.5 and 25°C.

Reagents: Phosphate-citrate-buffer 100 mM, pH 4.5
 ABTS 6 mM
 Enzyme solution ca. 0.1 U/ml
 All reagents should be prepared from deionized, distilled water. Prepare an enzyme solution in deionized, distilled water containing approximately 0.1 U/ml.

Procedure: *Pipette into a suitable quartz cuvette:*

buffer	700 μl
Enzyme solution	100 μl

Incubate at 25°C and check the temperature, place the reaction mixture into the photometer and monitor the absorbance until constant, start the enzymatic reaction by adding

ABTS-solution	200 μl
---------------	--------

Immediately mix and record the increase in absorbance at 420 nm for ca. 1 minute. Calculate the ΔA/min from the linear range of the curve (initial range: standardized as 10 sec after start).

Calculation:

$$\begin{aligned} \text{Volume activity} &= \frac{\Delta A/\text{min} \cdot 1000 \cdot \text{volume}_{\text{reaction mixture}} (\mu\text{l})}{36,000 \cdot \text{volume}_{\text{enzyme solution}} (\mu\text{l})} \\ &= [\text{U/ml enzyme solution}] \end{aligned}$$

References: [1] Childs RE & Bardsley WG (1975). *Biochem. J.* 145: 93-103.
 [2] Wolfenden BS & Willson RL (1982). *J. Chem. Soc. Perkin. Trans. 7:* 805-812.

Enzymatic Assay of Lignin-Peroxidase

Principle: Lignin peroxidase catalyzes the oxidation of non-phenolic aromatic rings into aryl cation radicals by H₂O₂. Aryl cation radicals are unstable and undergo various following reactions. A typical example is the oxidation of veratryl alcohol into veratryl aldehyde via the intermediary formation of veratryl cation and benzyl radicals [1].

Reaction: 3,4-dimethoxybenzyl alcohol + H₂O₂ → 3,4-dimethoxybenzaldehyde + 2 H₂O

Method: Continuous spectrophotometric determination of the increase in absorbance at 310 nm resulting from the formation of veratryl aldehyde ($\epsilon_{310\text{nm}} = 9,300 \text{ M}^{-1} \text{ cm}^{-1}$) [2].

Conditions: Temperature: 25°C
pH: 3.3
Wavelength: 310 nm
Light path: 1 cm

Unit definition: The amount of enzyme that oxidizes 1.0 µmole of 3,4-dimethoxybenzyl alcohol (veratryl alcohol) per minute at pH 3.3 at 25°C.

Reagents: Sodium-tartrate buffer 100 mM, pH 3.3
H₂O₂ 20 mM
Veratryl alcohol 50 mM
Enzyme solution ca. 0.5 U/ml

All reagents should be prepared from deionized, distilled water. Dissolve the enzyme at a concentration of 0.4 mg/ml in distilled water. Dilute further in distilled water if necessary.

Procedure: *Pipette into a suitable quartz cuvette:*

buffer	675 µl
veratryl alcohol	75 µl
Enzyme solution	200 µl

Incubate at 25°C and check the temperature, place the reaction mixture into the photometer and monitor the absorbance until constant, start the enzymatic reaction by adding

H ₂ O ₂	50 µl
-------------------------------	-------

Immediately mix and record the increase in absorbance at 270 nm for ca. 1 minute. Calculate the $\Delta A/\text{min}$ from the linear range of the curve (initial range: standardized as 10 sec after start).

Calculation:

$$\begin{aligned} \text{Volume activity} &= \frac{\Delta A/\text{min} \cdot 1000 \cdot \text{volume}_{\text{reaction mixture}} (\mu\text{l})}{9,300 \cdot \text{volume}_{\text{enzyme solution}} (\mu\text{l})} \\ &= [\text{U/ml enzyme solution}] \end{aligned}$$

References: [1] Schoemaker H et al. (1994). *FEMS Microbiol. Rev.* 13: 321-332.
[2] Tien M & Kirk TK (1984). *Proc. Nat. Acad. Sci. U.S.A.* 81: 2280-2284